

Luzerne Optical Laboratories Ltd

180 N. Wilkes-Barre Boulevard
Wilkes-Barre, PA 18702-5341
570-822-3183

Vision@LuzerneOptical.com
www.LuzerneOptical.com

Varilux Physio Enhanced™ Lenses **For the sharpest vision in any light, even at night**

Introducing Varilux Physio Enhanced with New W.A.V.E Technology 2™

- Superior aberration control
- Customized for optimal performance
- Digitally surfaced

THE NEXT GENERATION IN ABERRATION CONTROL

No other progressive addition lens comes close to minimizing distortion like the new Varilux Physio Enhanced lenses.

W.A.V.E. Technology 2™ is Essilor's unique next generation of lens design. It identifies and manages distortions other manufacturers miss-for the sharpest vision at every distance in every lighting condition, even in low light.

CUSTOMIZED FOR OPTIMAL PERFORMANCE

In addition to minimized distortion, Varilux Physio Enhanced lenses achieve a dramatically higher level of sharpness by customizing the design to three critical factors:

- Patient prescription
- Viewing distance-from near to far
- Light conditions-from bright light to low light

To ensure the most precise transfer of this new design to the lens, all Varilux Physio Enhanced lenses are digitally surfaced.

SHARPER VISION AT EVERY DISTANCE IN ANY LIGHT

Even at night, Varilux Physio Enhanced™ lenses preserve 93% of image sharpness.

ADVANCE YOUR PRACTICE WITH VARILUX PHYSIO ENHANCED

The sharper vision and reduced distortion enabled by W.A.V.E. Technology 2™ helps you provide practical patient benefits.

- Safer, more confident night driving
- Less strain during low-light reading, television and computer usage
- Reduction of swim and ease of adaptation.

Varilux Physio Enhanced lenses can strengthen your practice by helping you stand out from retail competition and decreasing non adapts and returns.

AVAILABILITY:

Available fitting height from 14 to 25+ mm

Product	Color	Power Range	Cylinder	Add Power
1.50 Index Plastic		+6.00D to -10.00D	-6.00	0.75-3.50,4.00
1.50 Transitions VI	Gray/Brown	+6.00D to -10.00D	-6.00	0.75-3.50,4.00
1.50 Xperio Polarized	Gray/Brown	+6.00D to -10.00D	-6.00	0.75-3.50
Airwear		+6.00D to -10.00D	-6.00	0.75-3.50,4.00
Airwear Transitions	Gray/Brown	+6.00D to -10.00D	-6.00	0.75-3.50,4.00
Airwear Xperio Polarized	Gray/Brown	+6.00D to -10.00D	-6.00	0.75-3.50
Thin & Lite 1.60		+8.00D to -12.00D	-6.00	0.75-3.50,4.00
Thin & Lite 1.60 Transitions VI	Gray/Brown	+8.00D to -12.00D	-6.00	0.75-3.50,4.00
Thin & Lite 1.67		+9.00D to -12.00D	-6.00	0.75-3.50,4.00
Thin & Lite 1.67 Transitions VI	Gray/Brown	+9.00D to -12.00D	-6.00	0.75-3.50,4.00
Thin & Lite 1.74*		+13.00D to -18.00D	-6.00	0.75-3.50,4.00

*Available only with Crizal Alize or Crizal Avance with Scotchgard Protector.

Essilor recommends Crizal Avance with Scotchgard Protector anti-reflective treatment on Varilux Physio Enhanced lenses to protect against glare, scratched, dust and fingerprints.