

EFFECTIVE | AUG 2010 | RELEASE.01

Eyecode[™] lenses take individualized vision to a new level by incorporating a patient's eye, frame, and behavior data. Patients will benefit from having the most precise vision, instantly and effortlessly.

VARILUX PHYSIO ENHANCED EYECODE[™]

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions [®] VI (Gray / Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio [®] (Gray / Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex [®]	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex [®] Transitions [®] VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear [®]	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear [®] Transitions [®] VI (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear [®] Transitions [®] VI XTRActive [™] (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear [®] Xperio [®] (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite [®] 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite [®] 1.60 Transitions [®] VI (Gray Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite [®] 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite [®] 1.67 Transitions [®] VI (Gray Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00

1 Minimum fitting height is 14mm. Note: At fitting heights less than 17mm, add power is up to +3.00 D.
Note: 1.74 will be available in September.

ESSILOR EYECODE[™] SINGLE VISION

Material	Power Range	Base Curves	Add Power	Cylinder
Airwear [®]	+8 to -10	N/A	N/A	-4.00
Airwear [®] Transitions [®] VI (Gray Brown)	+8 to -10	N/A	N/A	-4.00
Thin&Lite [®] 1.67	+8 to -13	N/A	N/A	-4.00
Thin&Lite [®] 1.67 Transitions [®] VI (Gray Brown)	+8 to -13	N/A	N/A	-4.00

Note: 1.74 will be available in September.

EFFECTIVE | AUG 2010 | RELEASE.01

Personalized lenses for superior performance. Changing conditions of wear create aberrations. Lens personalization assures optimal vision for each patient in each lens design.

VARILUX PHYSIO ENHANCED FIT™

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions® VI (Gray Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio® (Gray Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex®	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex® Transitions® VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear®	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI XTRActive™ (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Xperio® (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite® 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.60 Transitions® VI (Gray Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67 Transitions® VI (Gray Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.74 ¹	+13 to -18	1.75, 2.75, 3.75, 4.75, 6.00, 7.50	0.75 to 3.50, 4.00	-6.00

¹ Minimum fitting height is 14mm. Note: At fitting heights less than 17mm, add power is up to +3.00 D.
All lenses with Thin&Lite 1.74 comes with a choice of Crizal Sapphire™ with Scotchgard™ Protector, Crizal Avancé™ with Scotchgard™ Protector, or New and Improved Crizal® Alize®.

ESSILOR FIT™ SINGLE VISION

Material	Power Range	Base Curves	Add Power	Cylinder
Airwear®	+8 to -10	N/A	N/A	-4.00
Airwear® Transitions® VI (Gray Brown)	+8 to -10	N/A	N/A	-4.00
Thin&Lite® 1.67	+8 to -13	N/A	N/A	-4.00
Thin&Lite® 1.67 Transitions® VI (Gray Brown)	+8 to -13	N/A	N/A	-4.00
Thin&Lite® 1.74 ¹	+10 to -14	N/A	N/A	-4.00

¹ Minimum fitting height is 14mm. Note: At fitting heights less than 17mm, add power is up to +3.00 D.
All lenses with Thin&Lite 1.74 comes with a choice of Crizal Sapphire™ with Scotchgard™ Protector, Crizal Avancé™ with Scotchgard™ Protector, or New and Improved Crizal® Alize®.

For more information contact Luzerne Optical Laboratories LTD at 800-233-9637 or on the web at www.LuzerneOptical.com

EFFECTIVE | AUG 2010 | RELEASE.01

This breakthrough lens technology allows you to address the specific visual needs of Asian lens wearers in these key areas—eye anatomy, facial anatomy, and reading behavior.

VARILUX PHYSIO ENHANCED AZIO™

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions® VI (Gray Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio® (Gray Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex®	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex® Transitions® VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear®	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI XTRActive™ (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Xperio® (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite® 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.60 Transitions® VI (Gray Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67 Transitions® VI (Gray Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00

Note: 1.74 will be available in September.

ESSILOR AZIO™ SINGLE VISION

Material	Power Range	Base Curves	Add Power	Cylinder
Airwear®	+8 to -10	N/A	N/A	-4.00
Airwear® Transitions® VI (Gray Brown)	+8 to -10	N/A	N/A	-4.00
Thin&Lite® 1.67	+8 to -13	N/A	N/A	-4.00
Thin&Lite® 1.67 Transitions® VI (Gray Brown)	+8 to -13	N/A	N/A	-4.00

Note: 1.74 will be available in September.

EFFECTIVE | AUG 2010 | RELEASE.01

This breakthrough lens technology allows you to address the specific visual needs of Indian lens wearers in these key areas—eye anatomy, facial anatomy, and reading behavior.

VARILUX PHYSIO ENHANCED INDIA™

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions® VI (Gray Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio® (Gray Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex®	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex® Transitions® VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear®	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI XTRActive™ (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Xperio® (Gray Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite® 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.60 Transitions® VI (Gray Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67 Transitions® VI (Gray Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00

Note: 1.74 will be available in September.

For more information contact Luzerne Optical Laboratories LTD at 800-233-9637 or on the web at www.LuzerneOptical.com