

EFFECTIVE | AUG 2010 | RELEASE.01

Eyecode<sup>™</sup> lenses take individualized vision to a new level by incorporating a patient's eye, frame, and behavior data. Patients will benefit from having the most precise vision, instantly and effortlessly.

## VARILUX PHYSIO ENHANCED EYECODE<sup>™</sup>

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions <sup>®</sup> VI (Gray / Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio <sup>®</sup> (Gray / Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex <sup>®</sup>	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex <sup>®</sup> Transitions <sup>®</sup> VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear <sup>®</sup>	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear <sup>®</sup> Transitions <sup>®</sup> VI (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear <sup>®</sup> Transitions <sup>®</sup> VI XTRActive <sup>™</sup> (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear <sup>®</sup> Xperio <sup>®</sup> (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite <sup>®</sup> 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite <sup>®</sup> 1.60 Transitions <sup>®</sup> VI (Gray   Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite <sup>®</sup> 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite <sup>®</sup> 1.67 Transitions <sup>®</sup> VI (Gray   Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00

1 Minimum fitting height is 14mm. Note: At fitting heights less than 17mm, add power is up to +3.00 D.  
Note: 1.74 will be available in September.

## ESSILOR EYECODE<sup>™</sup> SINGLE VISION

Material	Power Range	Base Curves	Add Power	Cylinder
Airwear <sup>®</sup>	+8 to -10	N/A	N/A	-4.00
Airwear <sup>®</sup> Transitions <sup>®</sup> VI (Gray   Brown)	+8 to -10	N/A	N/A	-4.00
Thin&Lite <sup>®</sup> 1.67	+8 to -13	N/A	N/A	-4.00
Thin&Lite <sup>®</sup> 1.67 Transitions <sup>®</sup> VI (Gray   Brown)	+8 to -13	N/A	N/A	-4.00

Note: 1.74 will be available in September.

EFFECTIVE | AUG 2010 | RELEASE.01

Personalized lenses for superior performance. Changing conditions of wear create aberrations. Lens personalization assures optimal vision for each patient in each lens design.

## VARILUX PHYSIO ENHANCED FIT™

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions® VI (Gray   Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio® (Gray   Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex®	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex® Transitions® VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear®	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI XTRActive™ (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Xperio® (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite® 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.60 Transitions® VI (Gray   Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67 Transitions® VI (Gray   Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.74 <sup>1</sup>	+13 to -18	1.75, 2.75, 3.75, 4.75, 6.00, 7.50	0.75 to 3.50, 4.00	-6.00

<sup>1</sup> Minimum fitting height is 14mm. Note: At fitting heights less than 17mm, add power is up to +3.00 D.  
All lenses with Thin&Lite 1.74 comes with a choice of Crizal Sapphire™ with Scotchgard™ Protector, Crizal Avancé™ with Scotchgard™ Protector, or New and Improved Crizal® Alize®.

## ESSILOR FIT™ SINGLE VISION

Material	Power Range	Base Curves	Add Power	Cylinder
Airwear®	+8 to -10	N/A	N/A	-4.00
Airwear® Transitions® VI (Gray   Brown)	+8 to -10	N/A	N/A	-4.00
Thin&Lite® 1.67	+8 to -13	N/A	N/A	-4.00
Thin&Lite® 1.67 Transitions® VI (Gray   Brown)	+8 to -13	N/A	N/A	-4.00
Thin&Lite® 1.74 <sup>1</sup>	+10 to -14	N/A	N/A	-4.00

<sup>1</sup> Minimum fitting height is 14mm. Note: At fitting heights less than 17mm, add power is up to +3.00 D.  
All lenses with Thin&Lite 1.74 comes with a choice of Crizal Sapphire™ with Scotchgard™ Protector, Crizal Avancé™ with Scotchgard™ Protector, or New and Improved Crizal® Alize®.

For more information contact Luzerne Optical Laboratories LTD at 800-233-9637 or on the web at [www.LuzerneOptical.com](http://www.LuzerneOptical.com)

EFFECTIVE | AUG 2010 | RELEASE.01

This breakthrough lens technology allows you to address the specific visual needs of Asian lens wearers in these key areas—eye anatomy, facial anatomy, and reading behavior.

### VARILUX PHYSIO ENHANCED AZIO™

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions® VI (Gray   Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio® (Gray   Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex®	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex® Transitions® VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear®	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI XTRActive™ (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Xperio® (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite® 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.60 Transitions® VI (Gray   Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67 Transitions® VI (Gray   Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00

Note: 1.74 will be available in September.

### ESSILOR AZIO™ SINGLE VISION

Material	Power Range	Base Curves	Add Power	Cylinder
Airwear®	+8 to -10	N/A	N/A	-4.00
Airwear® Transitions® VI (Gray   Brown)	+8 to -10	N/A	N/A	-4.00
Thin&Lite® 1.67	+8 to -13	N/A	N/A	-4.00
Thin&Lite® 1.67 Transitions® VI (Gray   Brown)	+8 to -13	N/A	N/A	-4.00

Note: 1.74 will be available in September.

EFFECTIVE | AUG 2010 | RELEASE.01

This breakthrough lens technology allows you to address the specific visual needs of Indian lens wearers in these key areas—eye anatomy, facial anatomy, and reading behavior.

## VARILUX PHYSIO ENHANCED INDIA™

Material	Power Range	Base Curves	Add Power	Cylinder
1.50 Index Plastic	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Transitions® VI (Gray   Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
1.50 Xperio® (Gray   Brown)	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50	-6.00
Trivex®	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Trivex® Transitions® VI Gray	+6 to -10	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear®	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Transitions® VI XTRActive™ (Gray)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00	-6.00
Airwear® Xperio® (Gray   Brown)	+6 to -10	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00
Thin&Lite® 1.60	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.60 Transitions® VI (Gray   Brown)	+8 to -12	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00
Thin&Lite® 1.67 Transitions® VI (Gray   Brown)	+9 to -12	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00	-6.00

Note: 1.74 will be available in September.

For more information contact Luzerne Optical Laboratories LTD at 800-233-9637 or on the web at [www.LuzerneOptical.com](http://www.LuzerneOptical.com)